# ANIOS SPS PREMIUM

# ANIOS SPS Premium

# Sommaire

---

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### *Regulation and general information*

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*ANIOS SPS *Premium* is designed, produced and controlled by the Laboratoires ANIOS, certified by the AFAQ Organism under the number 1995/3723 in accordance with the ISO 9001 Quality System*

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*ANIOS SPS *Premium* is in compliance with the requirements of the European regulation related to biocides for disinfection of surfaces (Group 1 - PT 2)*

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*ANIOS SPS *Premium* formulation fulfils the requirements of Regulation (EC) No 648/2004 on detergents and its amendments*

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*ANIOS SPS *Premium* is labelled in accordance with the European regulation related to the classification and labelling of chemical products*

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*ANIOS SPS *Premium* meets our voluntary eco-design commitment included in our ANIOSAFE charter*

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### Pure product ready to use :

- *Storage between +5°C and +35°C.*
- *Stability:*
  - *Not opened spray: 3 years from the manufacturing date indicated on the label.*
  - *Opened spray: preserving of the 3 years if the product is stored in a container closed by the sprayer.*

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*Not concerned by biodegradability notion*

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*Readily biodegradable substances and/or inherently biodegradable*<sup>3</sup>

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*Substances which not fulfil the biodegradability criteria, or substances without any biodegradability data.*

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**ANIOS SPS PREMIUM** contains at least 99.7 % of inorganic material and of biodegradable organic material.

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##### Environment

*This mixture does not present an environmental hazard. No known or foreseeable environmental damage under standard conditions of use.\**

\* Classification drawn up in accordance with European regulation related to the classification and labelling of chemical products.

#### PERSONAL PROTECTION UNDER THE CONDITIONS OF THE INTENDED USE

*Refers exclusively to this product. If the product is handled and/or exposed simultaneously with other chemical agents, these must be taken into consideration when choosing personal safety equipment.*

##### Respiratory protection

*In case of insufficient ventilation, with the risk of exceeding the TLVs/TWA\*\*, wear suitable breathing apparatus.*

*Particularly breathing apparatus, A2P2 type.*

##### Hand protection

*When handling this product, wear suitable gloves.*

*Nitrile, latex or vinyl gloves.*

*Gloves must be replaced immediately if they show signs of wear and tear.*

##### Eye and face protection

*Avoid contact with eyes.*

*Before handling, wear safety goggles with protective sides accordance with standard EN166.*

*Prescription glasses are not considered as protection.*

*Provide eye washes in workshops.*

*If necessary, the nearest water supply*

\*\* please refer to SDS §8 for TLVs (Threshold Limit Values)/TWA (Time-Weighted Average) of concerned substances present in the product

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*The Safety Data Sheet (SDS) of this product is available on our Web site [www.anios.com](http://www.anios.com). In addition to visualization and/or impression, it is possible for you to download it.*

*The consultation of this SDS via internet enables you to receive later on by email any updated version of this document.*

*In the absence of an Internet access, the SDS remains available on request near our laboratories (tel.: +33 3 20 67 67 67 - fax: +33 3 20 67 67 68)*

\* an access is also possible at the site <http://www.quickfds.com>